

# Mobile Application

## Pairing Badge Shift User for iOS

### Pairing an Access Badge

1. Long press the PerfectServe Clinical Collaboration logo in the top left and select **Set Login Option**.
2. Select the appropriate secondary login option.
  - a. Barcode
  - b. NFC
3. Proceed to log into the PerfectServe Clinical Collaboration mobile app, select the three-line menu at the top left.
4. Select **Pair Access Card**.
5. When the camera opens, scan the barcode/NFC.
6. Hold the camera over until it is recognized.
7. Enter a **4-digit pin** and confirm your pin.
  - a. This is the pin you will enter every time you scan your badge to login.  
*Note: If you receive an error, please contact your Administration.*
8. You will receive a successful message window if your access badge has been added.
9. The next time you login, select **Tap here to scan barcode** for an easy badge access login.

